

Optical Cleaning Supplies

- Kimwipes®
- Cotton Swabs
- Pump Bottles
- Squeeze Bottles
- Finger Cots
- Latex Gloves



OptoSigma offers optical cleaning supplies to help you maintain your optics in pristine condition. Kimwipes and cotton swabs are perfect for use with alcohol or acetone to clean most delicate optical surfaces. We offer two types of refillable solvent bottles; pump and squeeze. Both are convenient for cleaning optics and are shipped empty. The choice is simply a matter of your preference. Latex finger cots and gloves will ensure that finger oils do not contaminate your optics during the cleaning process.

The following optical cleaning supplies will be found useful as accessories to your optical laboratory or production area.

Optical Cleaning Supplies

Description	Size	Quantity	Price	PART NUMBER
Kimwipes®	280 Sheets	1		098-0151
Cotton Swabs	100 Swabs	1		098-0161
Pump Bottle	8 oz	1		098-0141
Squeeze Bottle (Acetone)	8 oz	5		098-0131
Squeeze Bottle (Isopropanol)	8 oz	5		098-0132
Squeeze Bottle (Methanol)	8 oz	5		098-0133
Finger Cots	Small	1440		098-0111
Finger Cots	Medium	1440		098-0112
Finger Cots	Large	1440		098-0113
Latex Gloves	Small	100		098-0121
Latex Gloves	Medium	100		098-0122
Latex Gloves	Large	100		098-0123

ORDERING
&
TECHNICAL SUPPORT
(949) 851-5881
FAX (949) 851-5058
E-MAIL
sales@optosigma.com
WEB
www.optosigma.com

