

Mirror Mounts

- Optic Holders
- Optical Baseplates
- Optical Rails and Carriers
- Bases
- Rod Mounting Systems
- Post, Post Holders, Clamps

Rod Mounted Mirror Mounts

- Damped stainless steel rods for greater stability
- Compatible with existing 1.5 inch systems
- Precision adjustable mounts

These 1.5 inch diameter damped steel rods provide a very stable and vibration free method of mounting mirrors. The tilt mount is made to accept a range of replaceable OneTouch™ mirror holders.



Specifications & Tolerances

Dimensions: ±0.2mm
 Material: Rods, Stainless Steel
 Mounts, Aluminum

Rod Mounted Mirror Mounts

Item	Length, L		Price	INCH	METRIC
	(inches)	(mm)		PART NUMBER	PART NUMBER
Damped Rod	7.0	177.8		148-3310	148-3315
Damped Rod	14.0	355.6		148-3320	148-3325
Tilt Mount	-	-		112-5630	112-5630

OneTouch™ Mirror Adapters for Tilt Mount

Mirror Diameter	Price	PART NUMBER
(inches)		(mm)
1.0	25.4	112-5631
1.5	38.1	112-5632
2.0	50.8	112-5633

ORDERING
&
TECHNICAL SUPPORT
(949) 851-5881
FAX (949) 851-5058
E-MAIL
sales@optosigma.com
WEB
www.optosigma.com

