

Fixed Lens Holders for Multi-Elements

INCH **METRIC**

Diameter, D (mm)	Thread (mm)	Center Height, h (mm)	Price	PART NUMBER	PART NUMBER
56.0	M29 P0.75	27.0		111-0240	111-0245
60.0	M34 P0.75	29.0		111-0250	111-0255
70.0	M50.9 P0.75	34.0		111-0260	111-0265
56.0	M34.85 P0.75	27.0		111-0270	111-0275
66.0	M44.85 P0.75	32.0		111-0280	111-0285
76.0	M54.85 P0.75	37.0		111-0290	111-0295

Fixed Lens Holders - Part Number Compatibility Chart

INCH		METRIC		To Hold Multi-Element Lens	
PART NUMBER	PART NUMBER				
111-0240	111-0245	027-0210			
111-0250	111-0255	027-1210	027-0220	026-5430	
		027-1220	027-0230	026-5440	
		027-1225	027-0240	026-5442	
		027-1230	027-0250	026-5445	
		027-1240	027-0260	026-5450	
		027-1250	027-0270	026-5455	
		027-1260	027-0280	026-5460	
111-0260	111-0265	027-1350			
		027-1360			
		027-1370			
		027-1380			
		027-1390			
111-0270	111-0275	027-3010			
		027-3015			
		027-3020			
		027-3025			
		027-3030			
111-0280	111-0285	027-3035			
		027-3040			
		027-3045			
		027-3050			
		027-3055			
111-0290	111-0295	027-3060			
		027-3065			
		027-3070			
		027-3075			
		027-3080			

ORDERING
&
TECHNICAL SUPPORT
(949) 851-5881
FAX (949) 851-5058
E-MAIL
sales@optosigma.com
WEB
www.optosigma.com

