

Cross Clamps

- Half inch/half inch and half inch/6mm
- Perpendicular and rotational
- Strong clamping

Cross clamps are used to connect two posts together. The half inch posts may be joined to each other or to 6mm posts from the mini-system found on page 342. Both fixed perpendicular and rotational clamps are offered. These clamps are of the split ring type which provides a high clamping force. They are equipped with capstan screws for ease of use.



Cross Clamps

Type	Post Diameter (mm)	Price	PART NUMBER
Fixed	12.7 x 12.7		148-2350
Fixed	12.7 x 6.0		148-2340
Rotary	12.7 x 12.7		148-2250
Rotary	12.7 x 6.0		148-2240

ORDERING
&
TECHNICAL SUPPORT
(949) 851-5881
FAX (949) 851-5058
E-MAIL
sales@optosigma.com
WEB
www.optosigma.com

