

Metric Cross Clamps

- Perpendicular and rotational
- 12mm and 6mm, 12mm and 12mm,
- 20mm and 12mm, 20mm and 20mm
- Strong clamping

Cross clamps are used to connect two posts together of same or different diameters. Both fixed perpendicular and rotational clamps are offered. These clamps are of the split ring type which provides a high clamping force. They are equipped with capstan screws for ease of use.



Cross Clamps

Type	Post Diameter (mm)	Price	PART NUMBER
Fixed	12 x 6		148-2335
Fixed	12 x 12		148-2345
Fixed	20 x 12		148-2355
Fixed	20 x 20		148-2365
Rotary	12 x 6		148-2235
Rotary	12 x 12		148-2245
Rotary	20 x 12		148-2255
Rotary	20 x 20		148-2265