



Mini Translation Motorized Stages - 5 phase stepping motor

- Patented extended contact bearing rail with hardened steel body
- Two-plate, low-profile, high stiffness design
- Compatible with our existing five phase stepping motor stage controllers
- Remote control, manual joy stick control and micro stepping functions available



TSDM40-15X
Attached Screw
M3X8mmX4



TSDM60-20X
Attached Screw
M4X6mmX4

Specifications & Tolerances

Positional Repeatability: 0.010mm
 Backlash: 0.020mm
 Cables: MINI for PAT-001, SHOT-102
 DMINIS for SHOT-602, SHOT-202, SHOT-204MS
 Travel per pulse: 1 micron/0.5micron (Full/Half)

ORDERING
&
TECHNICAL SUPPORT

(949) 851-5881
FAX (949) 851-5058

E-MAIL
sales@optosigma.com

WEB
www.optosigma.com

Travel (mm)	Load Cpty. (Kg)	Axes of Travel	Weight (Kg)	Price	Model	PART NUMBER
15	3	X	0.35		TSDM40-15X	507-0040
20	4	X	0.6		TSDM60-20X	507-0060

