


Accessories

OptoSigma accessories include the plug and play control pad to be used with the SHOT-202 and SHOT-204MS to manually program and control the movement of the SGSP Series motorized stages.



The grease gun and grease tubes are for lubrication of the stages by the user at the appropriate time on site without returning to the factory. When used as recommended this extends the use of the product over a longer lifetime.



After selecting the stage and controller that meets your need, you will then select the required cables to meet your specific application. We offer a complete line of cables of various lengths and will customize to your requirements.



- Piezo Stages
- Linear Stages
- Rotation Stages
- ZF Stages
- Motorized Optics Holders
- Goniometers
- GS Closed Loop Stages
- Controllers & Drivers
- Cables & Accessories
- Software

Control Pads and Lubrication Gun

Description	Price	Model	PART NUMBER
Power adapter for PAT-001		PAT-001-POW1	555-5035-S
Control Pad for SHOT-202,204MS,SHOT-302GS,SHOT-304GS		CJ-200	555-4020
Joystick controller for SHOT-302GS,SHOT-304GS		JS-300	555-4021
Jog Operation box for SHOT-302GS,SHOT-304GS		JB-400	555-4022
Control Pad for SHOT-602		SCT-602	555-4040-S
Grease Gun		GREASE GUN	591-1110
Grease For SGSP20/SGSP26/Rotation Stages		AFA GREASE	591-1121
Grease For SGSP33/SGSP46/SGSP55 Stages		AFB GREASE	591-1122